

**BROEN-LAB** 

Table of Content / Product Overview	2
Introduction / Customer Statements	3

### **Emergency Showers & Eye Washes**

Hand Held Eye Showers & Accessories NEW	4-5
Eye Showers	6-7
Recessed Emergency Showers	8-9
Body Showers for Wall & Ceiling	10
Body & Eye Showers for Wall Mounting	11
Showers with Flexible Choice of Valve	11
Free Standing Body & Eye Showers	12
Accessories for Emergency Showers	13-14
Test Equipment for Emergency Showers / Eye Wash Bottles	15


### **Heated Emergency Showers & Eye Washes**



Introduction		16
Heated Combination Showers	NEW	17
Heated Body Showers	NEW	18
Heated Eye Showers	NEW	19
Mobile & Portable Emergency Showers	NEW	19
Heated Tank Showers	NEW	20
Accessories for Heated Showers	NEW	21
General Information		22
International Standards	2	3-24
Technical Information		25
Notes	2	6-27

### **PRODUCT OVERVIEW**

- Please refer to the actual page for product versions



**Premium** 

Our extensive knowledge allows us to continually develop and optimize our products so that they meet our customer needs and specifications. Safety, functionality and quality are an integral part of product design and flexibility.

### **BROEN-LAB - A SAFE CHOICE**

- Quality since 1948
- Complete product range
- Compliance with international standards
- Proven functionality world-wide
- Spare parts available for a guaranteed 10-year period

### **CUSTOMER STATEMENTS**



### IN THE SERVICE INDUSTRIES

"For more than four years Falck has had an excellent cooperation with BROEN-LAB, Our customers must be able to rely on Falck to guarantee the best solutions and consequently we only provide service for emergency showers that comply with the requirements of Danish Standard 15154 Parts I and II. BROEN-LAB products are extremely easy to install

and they are also clearly marked and easily serviced. Quite simply, they provide a better and safer solution for all parties."

Michael Hansen, Product Manager



### IN INDUSTRY

Carlsberg Denmark

"At Carlsberg we do not accept a compromise with regard to quality. We have therefore replaced all our emergency shower systems by products from BROEN-LAB so that they comply with the EU norm EN 15154.

We want to look well after the people working for us and it is our distinct opinion that our workplace must be a place where you feel well and safe."

Coordinating Work Environment Representative



### IN THE LABORATORY



"For quite some time our company has been in charge of the maintenance and service of installations at Akademiska Hus at Medicinberget in Gothenburg. This includes a number of laboratories which again means a considerable number of emergency shower systems supplied by BROEN-LAB.

We check the emergency showers regularly and are very content with BROEN-LAB's reliable design and construction which make them easy to inspect and test."

Anders Bergvall, Project Manager VVS-Teknik Mölndal



### IN THE EDUCATIONAL SECTOR

"At DTU we teach more than 7,000 students within very different fields - from the metal industry to microbiology - and within each specific field unique safety aspects must be considered. Therefore it is essential for us that we have a cooperation partner who supplies a quality product and who can also advice and guide us on the specific risk profile of each classroom."

Bjarne Andersen, HVAC Manager at DTU Campus Service

### nne pharmaplan in the medical industry



"As an international advisor it is important for us to have a partner that complies with both national and international norms and standards - and who can document this by an independent third party certification. At the same it is important to have a partner that can offer customized solutions supporting the architectonic expression of the project instead of colliding with it. Even though we are talking about emergency equipment, beautiful design is a plus."

Henrik Als Nielsen. Engineer, Ventilation & Building Utilities, NNE Pharmaplan

BROEN-LAB has developed new hand-held eye showers, which comply to EN15154-2 and ANSI Z358.1 norms and also offer many improved features.



### TECHNICAL SPECIFICATIONS

Min. working pressure: 1.5 bar Max. working pressure: 10 bar Min. temperature: 5°C Max. temperature: 80°C

### MATERIALS WITH MEDIA CONTACT

Plastics: PA, POM, TPE Metals: SS, AISI 316 Rubber: EPDM

Materials not in media contact: PP, PA, SEBS

### SAFETY

- Improved trigger lock
- Improved activation of the shower



### **OPERATION**

- Stay-open dust cap
- Light weight
- Well-balanced construction
- Better spray pattern



### RELIABILITY

- Easy-to-clean spray head
- Anti-limescale design and material choice
- Improved functionality
- Longer lifetime





### Eye Shower with Two 45° Heads

### No. 27 10 009 10 05 011

Weight: Approx. 0.5 kg

Table mounted eye shower with two 45° angled heads, incl. 1.5m inlet hose. Trigger operated including a built-in check valve and FLOWFIX for regulation of water flow (14 l/min).





### Eye Shower with Two 45° Heads

### No. 27 10 009 10 03 011

Weight: Approx. 0.5 kg

Wall mounted eye shower with two 45° angled heads, incl. 1.5m inlet hose.
Trigger operated including a built-in check valve and FLOWFIX for regulation of water flow (14 I/min).





### Eye Shower with a Single 45° Head

### No. 27 10 003 10 05 011

Weight: Approx. 0.5 kg

Table mounted eye shower with a single 45° head, incl. 1.5m inlet hose. Trigger operated including a built-in check valve and FLOWFIX for regulation of water flow (7 l/min).





### Eye Shower with a Single 45° Head

### No. 27 10 003 10 03 011

Weight: Approx. 0.5 kg

Wall mounted eye shower with a single 45° head, incl. 1.5m inlet hose. Trigger operated including a built-in check valve and FLOWFIX for regulation of water flow (7 l/min).





Other configurations are possible - please refer to ordering information on page 5.

### **Shower for General Purpose Rinsing**

### No. 27 10 001 10 00 001

Weight: Approx. 0.1 kg

Shower for general purpose rinsing. Trigger operated including a built-in check valve. Without FLOWFIX for regulation.

Also available as straight model.

### Eye Shower with a Single 45° Head

### No. 27 10 003 10 00 001

Weight: Approx. 0.2 kg

Eye shower with a single head. Trigger operated including a built-in check valve and FLOWFIX for regulation of water flow (7 l/min).

Also available as straight model.

### Eye Shower with Two 45° Heads



### No. 27 10 009 10 00 001

Weight: Approx. 0.2 kg

Eye shower with two angled heads. Trigger operated including a built-in check valve and FLOWFIX for regulation of water flow (14 l/min).

Also available as straight model.



### **ACCESSORIES FOR HAND-HELD EYE SHOWERS**

### **Mounting Accessories (YY):**

### No. 17 760 300

Table base (M30x1.5, 60mm)



### No. 17 760 310

30° Wall base





### No. 17 760 320 Wall bracket

No. 17 760 330 Bracket for modular body shower



### Signs and Instructions (YY):

### No. 17 001 970

Complete set with signs and instructions

for eye shower



### No. 96M0339

Sign for eye shower (sticker)



### Hoses and hose accessories (ZZ):

### No. 93H40431500

1.5m Silicone hose with a Stainless Steel braiding for hand-held eye showers.



Length 2m: No. 93H40432000 (other lengths on request).

### No. 17 701 750

Counter weight/stop for hose



**BROEN-LAB** EMERGENCY SHOWER SYSTE

### ORDERING INFORMATION

For ordering a complete hand-held eye shower including accessories, please follow the following structure:

### 27 10 XXX 10 YY ZZ

### Serial number

- 001 Rinsing model, 45°
- 003 Eye Shower, single, 45°
- 004 Rinsing model, straight
- 006 Eye shower, single, straight
- 009 Eye shower, double, 45°
- 012 Eye shower, double, straight

### **Country code**

10 - International

### \* Accessories

- 00 None
- 01 17 001 970
- 02 96M0339
- 03 17 001 970 + Wall bracket
- 04 17 001 970 + 30° Wall Base
- 05 17 001 970 + Table base
- 06 96M0339 + Wall bracket
- 07 96M0339 + 30° Wall Base
- 08 96M0339 + Table base
- 09 96M0339 + Bracket for
  - modular body shower

### \* Hose length + Hose accessory

- 00 None
- 01 1.5m with G1/2 union nut
- 02 1.5m with G1/2 union nut + stop
- 03 2m with G1/2 union nut
- 04 2m with G1/2 union nut + stop

\* Please request document No. 99G0045 for a complete list of available hoses and accessories. For more information and technical datasheets please contact your BROEN-LAB sales representative.

The European Standard EN15154 advise that the water temperature delivered by emergency equipment should be "tepid".


Consequently, BROEN-LAB recommends the installation of a Thermostatic Mixing Valve (TMV) with safety shower equipment.

Please refer to page 13 for further information.



### Eye Shower, Wall Mounted





Weight: Approx. 1 kg

Eye shower, wall mounted with two heads in brass with chemical resistant BROEN-LAB Polycoat and with built-in BALLOFIX® ball valve for regulation of water amount.


Min. working pressure 1.5 bar

Min. working pressure: 1.5 bar.

### **Eye Shower, Wall Mounted**



### No. 17 069 009




### Weight: Approx. 2 kg

Eye shower, wall mounted for exposed pipe work in brass with chemical resistant BROEN-LAB Polycoat with built-in FLOWFIX for water regulation (26 l/min).

Min. working pressure: 1.5 bar.

### Eye/Face Shower with Bowl, Wall Mounted






### **Eye Shower, Wall Mounted**



No. 17 250 009



Weight: Approx. 4.5 kg


Eye shower, wall mounted in brass with chemical resistant BROEN-LAB Polycoat. The water is regulated with a built-in FLOWFIX (26 I/min).

Min. working pressure: 1.5 bar.

### Eye Shower, Table Mounted



No. 17 200 009






Weight: Approx. 2 kg

Eye shower, table mounted with two heads in brass with chemical resistant BROEN-LAB. Water is regulated with a built-in FLOWFIX (26 I/min). Min. working pressure: 1.5 bar.

### Eye Shower, Table Mounted, Swing Down






Weight: Approx. 4.5 kg

Eye shower, table mounted, swing down in brass with chemical resistant BROEN-LAB Polycoat for mounting on back side of the sink. The water amount is regulated with built-in FLOWFIX (26 l/min). Min. working pressure: 2 bar.

### Eye Shower, Table Mounted, Swivel Joint





Weight: Approx. 4 kg


Eye shower, table mounted, swivel joint in brass with chemical resistant BROEN-LAB Polycoat for mounting on the right side of sink. The water amount is regulated with the built-in FLOWFIX (26 l/min).

Min. working pressure: 1.5 bar.

## BROEN-LAB EMERGENCY SHOWER SYSTEM

### Eye/Face Shower with Bowl, Recessed





Eye/face shower with two heads in brass with chemical resistant BROEN-LAB Polycoat and Stainless Steel bowl for recessed mounting on table top, with push-down control. The water amount is regulated with built-in FLOWFIX (26 l/min). Min. working pressure: 1.5 bar.

### **SPECIAL FEATURES**

### **GENERAL:**

- Choose BROEN-LAB Polycoated surface or Stainless Steel.
- Cabinet and cabinet door in Stainless Steel AISI 304.
- · Disability-Friendly models comply with ANSI A117.1-2003.
- Easy to mount and install.
- · Cabinet can also be surface mounted (except 17 155 009).

Choose among different body shower outlets by ordering recessed units and outlets separately to fit your exact needs of installation and application.


E.g.: 17 154 009 + 17 110 009

### **BODY SHOWER:**

- Strong integrated handle for operating body shower.
- Free choice of combination between ceiling or wall mounted body shower outlet.
- · Easy to mount and install.
- Premounted valve for body shower.

### **Body Shower Outlet, Ceiling Mounted, Concealed**






Weight: Approx. 1 kg

Body shower outlet, ceiling mounted, concealed in Stainless Steel with selfdraining shower head. Inlet G 1" male. Min. working pressure: 1.5 bar.


**Body Shower Outlet, Wall Mounted** 



Weight: Approx. 1.5 kg

Body shower outlet, wall mounted, in brass with chemical resistant BROEN-LAB Polycoat with self-draining shower head. Inlet G 3/4" male. With adjustable flange.

Min. working pressure: 1.5 bar.



### **EYE WASH:**

- Cabinet and cabinet door in Stainless Steel AISI 304.
- Integrated drain (on models with door).
- Cabinet door catches the water and prevents it from splashing on the floor.
- Filter BALLOFIX® ensures clean shower heads.

### **Body Shower Outlet, Ceiling Mounted**






Body shower outlet, ceiling mounted, in brass with chemical resistant BROEN-LAB Polycoat with self-draining shower head. Fits different ceiling heights by adjusting the flange. Inlet G 3/4" male.

Min. working pressure: 1.5 bar.



Male inlet and fixed flange





Weight: Approx. 1.5 kg

Wall mounted body shower outlet in Stainless Steel AISI 304, with self-draining shower head. Inlet G 1" male. With adjustable flange. Dimensions comply with the requirements of ANSI A117.1-2003. Min. working pressure: 1.5 bar.

### Eye Shower, Wall Recessed



### **Eye Shower, Wall Recessed**





**Disability-Friendly** 

No. 17 152 009 Other Models


No 17 151 009


### Other Models

No. 17 150 009

### Anti-corrosive

Same as ill. Polycoat but anti-corrosive cabinet. No. 17 160 009





### Without cabinet door

### Anti-corrosive

Same as ill. Polycoat but anti-corrosive cabinet. No. 17 162 009

Without cabinet door No. 17 161 009

### Weight: Approx. 8 kg

Eye shower, wall recessed with cabinet with chemical resistant BROEN-LAB Polycoat, two outlet heads and built-in FLOWFIX (14 I/min).

Cabinet door activates the eye wash and catches the water.

Min. working pressure: 2 bar.

### Weight: Approx. 10 kg

Eye shower, wall recessed, disability-friendly with cabinet with chemical resistant BROEN-LAB Polycoat, two outlet heads and built-in FLOWFIX (14 I/min). Cabinet door activates the eye wash and catches the water. Min. working pressure: 2 bar.

### Eye and Body Shower, Wall Recessed







### **Disability-Friendly**

No. 17 154 009

### Other Models



### Anti-corrosive Same as ill. Polycoat but anti-corrosive cabinet.

Without cabinet door

No. 17 164 009

Water outlet DIA 1 1/2' (Ø40 mm) drain pipe

female inlet for body shower

### No. 17 163 009

### Weight: Approx. 16 kg


Eye shower, with 1 "valve for body shower, wall recessed, disability-friendly with cabinet with chemical resistant BROEN-LAB Polycoat, two outlet heads and built-in FLOWFIX (14 I/min). Cabinet door activates the eye wash and catches the water. Integrated hand activating 1" valve for body shower with various options for shower outlet. Min. working pressure: 2 bar.

### **Body Shower, Wall Recessed, Hand Activator**



### Other Models

### Anti-corrosive Same as ill. Polycoat but anti-corrosive cabinet. No. 17 165 009



### Weight: Approx. 2.5 kg

Wall recessed cabinet with chemical resistant BROEN-LAB Polycoat and hand activator for body shower with various options for shower outlet.

1" valve for 3/4" or 1" pipe.

Min. working pressure: 1 bar.

G 1/2"

G 1/2"

female

800 mm

(see page 13)

### **Body Shower, Ceiling Mounted** G 1/2" mal e

### No. 17 080 009 60 mm Weight: Approx. 2 kg 248 mm Body shower, ceiling mounted with self-draining head. In brass with chemical resistant

**BROFN-I AB** 

Min. working pres-

Polycoat.

sure: 1 bar.

600 mm

490 m

800 mm




### **Body Shower, Wall Mounted**

No. 17 973 009

Weight: Approx. 5.5 kg


Body shower, wall mounted with selfdraining shower head. In brass with chemical resistant BROEN-LAB Polycoat. Min. working pressure: 1 bar.



### **Body Shower, Wall Mounted**

No. 17 081 009

Weight: Approx. 3 kg Body shower, wall mounted with self-draining shower head. In brass with chemical resistant BROEN-LAB Polycoat. Min. working pressure: 1 bar.



Anti-corrosive - Modular

Body shower, wall mounted, modular in SS & Stanniloy with self-draining shower head for flexible wall mounting with pull rod.

Min. working pressure: 1 bar.



No. 17 970 004

Supporting wall flange can be used for mounting (see page 13)

660 mm

490 mm

### **Body Shower for Exposed Pipework, Wall Mounted**

### No. 17 082 009

Weight: Approx. 3 kg Body shower for exposed pipework with self-draining head. In brass with chemical resistant BROEN-LAB Polycoat. Min. working pressure: 1 bar.

### **Concealed Pipework**



Anti-corrosive - Modular

No. 17 971 004

**Disability-Friendly** female inlet G1" 584 mm 150 mm 1290 mm

Weight: 2.5 kg

Body shower, concealed, modular in SS with selfdraining shower head for recessed flexible wall mounting. With extended pull rod for disabled people. Min. working pressure: 1 bar.

### **Body Shower, Wall Mounted above Door**

Inlet female G1" (G3/4" female adapter incl.)

Supporting wall flange can be used for mounting

### No. 17 087 009

Weight: Approx. 4.5 kg Body shower, wall mounted above door for exposed pipework with self-draining head. In brass with chemical resistant BROEN-LAB Polycoat. Valve to be mounted left or right of door. Min. working pressure: 1 bar.

### female inlet G 1 500 mm 4 x ø10.5 mm 591 150 mm

### Anti-corrosive - Modular

No. 17 972 004


Weight: Approx. 6.5 kg Body shower, wall mounted above door for exposed pipework with self-draining head. In Stainless Steel and Stanniloy. Valve to be mounted left or right of door.

Min. working pressure: 1 bar.



500 mm


### **Body and Eye Shower, Wall Mounted**



chemical resistant BROEN-LAB Polycoat with self-draining shower head. Water inlet at the top or bottom. Eye shower has built-in FLOWFIX (26 I/min).

Min. working pressure: 1.5 bar.

### Body and Eye Shower, Wall Mounted, Modular



Body and eye shower, wall mounted, modular. In brass and SS with chemical resistant red BROEN-LAB Polycoat. Body shower has self-draining shower head. Water inlet at the top or bottom. Eye shower has built-in FLOWFIX (26 I/min). Min. working pressure: 1.5 bar.


### No. 17 952 009 Other Models



Same as ill. Polycoat but anti-corrosive. No. 17 963 004

With eye/face shower and bowl No. 17 964 004


### Body Shower, Wall mounted, Modular



Body shower, wall mounted, modular. In brass and SS with chemical resistant red BROEN-LAB Polycoat with self-draining shower head for flexible wall mounting with pull rod


Min. working pressure: 1 bar.

### **Body Shower, Concealed, Modular**



Body shower, concealed, modular. In SS and Chrome/Nickel plated brass with chemical resistant red BROEN-LAB Polycoat and self-draining shower head. For recessed flexible wall mounting. With extended pull rod for disabled people. Min. working pressure: 1 bar.

### SHOWERS WITH FLEXIBLE CHOICE OF VALVE



Body shower outlet for flexible choice of valve in brass and chemical resistant BROEN-LAB Polycoat.

Min. working pressure: 1 bar.

### Weight: Approx. 1 kg

85 mm

105 mr

Weight: Approx. 2.5 kg

BALLOFIX® ball valve


G 3/4" female

No. 17 022 009

190 mm

Manually hand operated of 1" BALLOFIX® ball valve for 3/4" pipe work.

### BALLOFIX® ball valve w./rod




Manually rod operated 1" BALLOFIX° ball valve for 3/4" pipe work.

Please refer to page 8 for more outlets for flexible choice of valve.

### Body Shower with Eye/Face Shower and Bowl, Modular


### No. 17 956 009 650 mm Other Models No. 17 955 009 inlet G1 With eve shower No. 17 954 009 nlet for Basic Model shower 2223 mm No. 17 956 609 With hand-held eye 1587 female outlet G1' shower 000 Weight: Approx. 17 kg Body shower with eye/face shower and bowl, modular In brass and SS with chemical resistant red BROEN-LAB Polycoat with self-draining shower head. Eye wash with built-in FLOWFIX for regulation of water flow (26 $\,$ I/min). Bowl in Stainless Steel with integrated outlet. Min. working pressure: 1.5 bar.


### Body Shower with Eye/Face Shower and Bowl, Modular



Body shower with eye/face shower and bowl, modular. In brass and SS with chemical resistant red BROEN-LAB Polycoat with self-draining shower head. Eye wash with built-in FLOWFIX for regulation of water flow (26 l/min). Bowl in Stainless Steel with integrated outlet.

Min. working pressure: 1.5 bar.






### Eye Shower, Free Standing, Modular






No. 17 950 009

No. 17 961 004



Weight: Approx. 7.5 kg
Eye shower, free standing,
modular in brass and SS
with chemical resistant red
BROEN-LAB Polycoat.
Built-in FLOWFIX for
regulation of water flow
(26 I/min). Min. working
pressure: 1.5 bar.





Other Models

No. 17 962 004
Eye/face shower
with bowl

female
inlet G1"

### Thermostatic Mixing Valves for Eye and Body Showers

BROEN-LAB TMVs for eye and body showers are thermostatic mixing valves designed for delivering tepid water (20-25°C) for eye/ hand or body showers.


The solution is designed for the operation of one eye/hand shower respectively one body shower.

The mixed temperature is adjusted to tepid water from factory but can easily be adjusted on site by installer or operator.

### The TMV incorporates in one body the following features:

- Fail-safe on hot side: If the cold water is interrupted the built-in anti-scalding device will automatically cease the flow of hot water. **AUTOMATIC OPERATION**
- Fail-safe on cold side: If the hot water is interrupted the built-in intentional leak will allow cold water to flow at a sufficient flowrate to ensure operation. THIS WORKS AUTOMATICALLY AS WELL
- · Including check valves on both hot and cold supply

### Thermostatic Mixing Valves for Eye Shower



Weight: Approx. 1.5 kg

### No. 17 140 009:

Thermostatic mixing valve for eye shower. Min. working pressure: 1 bar.

### No. 26 140 059:

Lockable cap for 17 140.009.

### **Thermostatic Mixing Valves for Body Shower**



### No. 17 459 009\*:

Thermostatic mixing valve for body shower. Min. working pressure: 1.5 bar.

No. 17 460 009\*: Model with G1" inlet/outlet connections

\*Please note: Reduced temperature stability at low flow rates!

### **Signs for Eye and Body Showers**

### Nr. 17 001 920

(Weight: Approx. 0.23 kg)

Sign (pull rod mounted) and instruction card for Body and Eye Shower

Sign (pull rod mounted) and

instruction card for Body Shower



### Nr. 17 002 27

(Weight: Approx. 0.03 kg) Instruction and Service Sign

**Body Shower** 



### Nr. 17 002 28

Nr. 17 002 29

(Weight: Approx. 0.03 kg)

Eye Shower fountain

(Weight: Approx. 0.03 kg) Instruction and Service Sign -

**Body and Eye Shower** 



(sticker)

Nr. 96M0042

Nr. 96M0043

(Weight: Approx. 0.01 kg)

Sign for body shower

(Weight: Approx. 0.01 kg) Sign for eye shower (sticker)



### Nr. 96M0345

(Weight: Approx. 0.01 kg) Sign for direction (sticker)



### Nr. 17 002 010

Nr. 17 001 910

(Weight: Approx. 0.13 kg)

(Weight: Approx. 0.15 kg)

Sign (wall mounted) and instruction card for Body and Eye Shower



### Nr. 17 002 30

(Weight: Approx. 0.03 kg) Instruction and Service Sign -

Instruction and Service Sign -



Eye Shower - Hand-Held



### Nr. 17 002 31

(Weight: Approx. 0.03 kg) Instruction and Service Sign -Body Shower with Hand-Held Eye Shower



### Wall bracket for body showers




### No. 17 001 079

Wall bracket for support of wall mounted body showers (size 1/2"). Other sizes on request

### Nr. 17 001 900

(Weight: Approx. 0.08 kg) Sign (wall mounted) and instruction card for Eye Shower




### Nr. 17 002 000 (Weight: Approx. 0.08 kg)

Sign (wall mounted) and instruction card for Body Shower





### **Cable**



No. 17 901 64

A cable may be installed as a supplement to the pull rod solution i.e. the cable enables the user to activate the shower at floor level.

### Foot Panel for Operation of Free Standing Body Shower



Weight: Approx. 23 kg

No. 17 919 009 - Stainless Steel No. 17 834 009 - Galvanized Steel

Foot panel for operation of body showers. Possible to combine with 17 909.909 if mounted 180° opposite.

Delivered with all mounting accessories

### **Alarm for Emergency Shower**



No. 17 851 95

Signal and light alarm for emergency showers.

- 112dB(A) output
- IP66 protection
- 32 alarm tones

Working temperature: -20°C to 70°C Power supply: 230VAC, 50/60Hz

For operation in Non Ex Proof areas.

### Flow switch



96P0062	brass	G3/4	25
96P0063	brass	G 1 1/4	25
96P0066	SS	G3/4	25
96P0068	SS	G 1 1/4	25

Item no. | Material | Size\*

Weight: Approx. 0,5 kg


\* Please note that size of the flow switch should be a size bigger than the inlet of the emergency shower

- Good reliability
- Low presure loss
- Dirt resistant
- Hermetic separation of electrical and hydraulic components
- Stress free fixation of switch unit by plastic cap

PN, bar

When there is a stream of water in the connection piping to the emergency shower the paddle in the flow switch is triggering an electrical reed switch. This switch can be connected to any type of warning system, such as an alarm central, a horn, a light or similar.

### Foot Treadle for Operation of Eye Showers



Weight: Approx. 1 kg

No. 17 909 909

Foot treadle in aluminium and SS for operation of eye showers. Possible to combine with 17 919.009 if mounted 180° opposite. Delivered with all mounting accessories.

### **Inlet hoses for Emergency Showers**



No. 17 701 72 (Male ISO 7/1-R1" - Female G1 union nut connection) No. 17 701 73 (Male ISO 7/1-R1" - OD28mm connection) No. 17 701 74 (Male G3/4 - Female G1 union nut connection)

Material: Silicone Braiding: Stainless Steel Inside diameter: Ø25mm Outside diameter: Ø34mm Max.working pressure: 10 bar Max.working temperature: 90°C


For more information and technical datasheets please contact your BROEN-LAB sales representative.

### **Test Trolley for Emergency Showers**

No. 17 851 92 - model with integrated pump No. 17 851 93 - model without pump

Test trolley with or without integrated pump. Drainage hose is included. Capacity: 90 I Can be used for heated and unheated emergency showers.

Weight: Approx. 22 kg



### **Test Funnel for Emergency Showers**

### No. 17 001 08

Test funnel for testing of emergency showers in PC (polycarbonate) with 2.2m hose.

### No. 27 99 008 10 01

Measuring Cylinder for checking Eye Wash Stations





### **EYE WASH BOTTLES**

### Duo Eye Wash Bottle, 500ml pH Neutral

### No. P7 17 003 10 00 000

DUO 500ml pH neutral.

Sterile 4.9 % Phosphate buffer solution.

For neutralizing acids and alkalis. Complete with an eye device for washing both eyes at the same time.

Durability 3 years.

Full box with 6 units available.



### Eye Wash Bottle (single eye), 200ml pH Neutral

### No. P7 17 004 10 00 000

200ml pH neutral.

A small, handy bottle containing 4.9% sterile phosphate buffer solution that neutralizes the eye's fluid after accidents involving acids or alkali.

The bottle is provided with an ergonomic eye cup, a dust cap and clear instructions on the label.

Eye Wash Bottle (single eye), 500ml Sterile Solution

Durability 3 years.

Full box with 10 units available.



### **Duo Eye Wash Bottle, 500ml Sterile Solution**

### No. P7 17 002 10 00 000

DUO 500ml Sterile.

Bottle containing a 0.9% sterile sodium chloride solution. Complete with an eye device for washing both eyes at the same time.

Durability 3 years.

Full box with 12 units available.



### No. P7 17 001 10 00 000

500ml Sterile.

Bottle containing a 0.9% sterile sodium chloride solution.

The bottle is provided with an ergonomic eye cup, a dust cap and clear instructions on the label.

Durability 3 years.

Full box with 12 units available.



### **Eye Wash Station**

### No. P7 33 001 10 00 402

Dust-proof wall mounted eye wash station in ABS plastic for two 500ml eye wash bottles + one 200ml eye wash bottle.

Especially suitable for dirty and dusty production areas and mobile work places.



### Eye Wash Station (classic model)

### No. P7 33 002 10 00 402

Dust-proof wall mounted eye wash station in polystyrene for two 500ml eye wash bottles. Especially suitable for dirty and dusty production areas and mobile work places.



### Wall Mount for Eye Washes (Duo or single eye)

### No. P7 53 001 10 00 409

Wall mount for pH neutral solutions (blue)

**No. P7 53 001 10 00 402**Wall mount for sterile

solutions (green)





### Wall Mount for Eye Washes (single eye)

### No. P7 53 002 10 00 409

Wall mount for pH neutral solutions (blue)

**No. P7 53 002 10 00 402** Wall mount for sterile solutions (green)





For more information and technical datasheets please contact your BROEN-LAB sales representative.

### INTRODUCTION

Thousands of industrial accidents happen every year in spite of precautionary safety measures, accidents that are caused by human exposure to heat (fire, boiling water and oil, steam, etc.), to corrosive agents (fluid acids and caustics) and/or to influence by dust and corrosive aerosols (metal, dirt, fluids).



These accidents cause the injured person acute pain, disablement and may even lead to loss of life. However, experience shows that the consequences and human suffering following such accidents can be considerably limited if water is applied directly to the affected areas instantly after the exposure.

It is a fact of life that the risk lurks in many places: in the laboratory, in the kitchen, in the warehouse, in the chemical plant, in the workshop or generally speaking in any place where heat, corrosion or dust can cause serious human injury. For these reasons we strongly recommend the installation of adequate emergency shower systems as a means of limiting these risks.

### Location and Operation of Emergency Showers

As a basic rule the shower must be located as close to the workplace as possible so that the general first aid instructions of flushing the affected area immediately after the accident has occurred can be followed.

Please note that if strong acids or caustics are involved immediate action is demanded. When installing the emergency shower system please also be aware that it must be easy to locate and access at all times for potential users.

BROEN-LAB Emergency Showers can be recognized by their red and white stripes on body and eye showers, red pull rods and luminescent signs and symbols for the clear identification of emergency shower



### Selecting the Right Shower

There are many different conditions that determine the type of emergency shower to be installed. It is important that the chosen unit is designed and tested according to the existing standard as well as suits the specific application and location - *and environment!* 



Tank showers - see page 21

In this section you can find a wide range of BROEN-LAB Emergency Showers that can be installed outside and be used in areas where outside temperatures go as low as -35°C.



Frost protected outdoor combination safety showers, trace tape heated and pre-insulated. Delivered with eye/facewash fountain and additional eyebath diffuser on flexible hose to treat unconsciuos or disabled personnel.


Available in Galvanized or Stainless Steel.

### **Heated Body Shower**

### - Free Standing with Eye Wash Fountain

## Self-draining PVC spray nozzle (60° angle) Eye/Facewash fountain Junction box

### Heated Body Shower - Free Standing with Eye Wash Fountain



### TECHNICAL SPECIFICATIONS

Min. working pressure: 2 Bar Min. working temperature: -20°C Water Flow: 115 I/min.

Supply voltage: 230 VAC, 50/60 Hz

Heated underground extension pipe (ordered separately, see page 22 for more information)

### TECHNICAL SPECIFICATIONS

Min. working pressure: 2 Bar Min. working temperature: -35°C  $\,$ 

Water Flow: 115 l/min.

Supply voltage: 230VAC, 50/60Hz

Heated underground extension pipe (ordered separately, see page 22 for more information)



### ORDERING INFORMATION

BROEN-LAB item no.	Material		Min. working temperature	Eye wash fountain	Mounting	Heated underground extension pipe
H7 41 002 10 00 321	Galvanized Steel	Yes	-20°C	Yes	Floor	Available
H7 41 002 10 00 221	Stainless Steel	Yes	-20°C	Yes	Floor	Available
H7 41 001 10 00 341	Galvanized Steel	Yes	-35°C	Yes	Floor	Available
H7 41 001 10 00 241	Stainless Steel	Yes	-35°C	Yes	Floor	Available
H7 41 001 10 00 331	Galvanized Steel	No	-35°C	Yes	Floor	Available
H7 41 001 10 00 231	Stainless Steel	No	-35°C	Yes	Floor	Available

<sup>\*</sup>Ex-proof zone 1 and zone 2.


For more information and technical datasheets please contact your BROEN-LAB sales representative.

Frost protected outdoor combination safety showers, trace tape heated and pre-insulated. Available in Galvanized or Stainless Steel.


### Heated Body Shower - Free Standing

# Self-draining PVC spray nozzle (60° angle) Shower pull rod Junction box Heated underground extension pipe (ordered separately, see page 22 for more information)

### Heated Body Shower - Free Standing



### Heated Body Shower - Wall Mounted



### TECHNICAL SPECIFICATIONS

Min. working pressure: 2 Bar Min. working temperature: -20°C

Water Flow: 115 l/min.

Supply voltage: 230VAC, 50/60Hz

### **TECHNICAL SPECIFICATIONS**

Min. working pressure: 2 Bar Min. working temperature: -35°C

Water Flow: 115 l/min.

Supply voltage: 230VAC, 50/60Hz

### TECHNICAL SPECIFICATIONS

Min. working pressure: 2 Bar Min. working temperature: -20°C

Water Flow: 115 I/min.

Supply voltage: 230VAC, 50/60Hz

### ORDERING INFORMATION

BROEN-LAB item no.	Material	Ex-proof*	Min. working temperature	Eye wash fountain	Mounting	Heated underground extension pipe
H7 21 002 10 00 321	Galvanized Steel	Yes	-20°C	No	Floor	Available
H7 21 002 10 00 221	Stainless Steel	Yes	-20°C	No	Floor	Available
H7 21 001 10 00 341	Galvanized Steel	Yes	-35°C	No	Floor	Available
H7 21 001 10 00 241	Stainless Steel	Yes	-35°C	No	Floor	Available
H7 21 001 10 00 331	Galvanized Steel	No	-35°C	No	Floor	Available
H7 21 001 10 00 231	Stainless Steel	No	-35°C	No	Floor	Available
H7 23 001 10 00 321	Galvanized Steel	Yes	-20°C	-	Wall	-
H7 23 001 10 00 221	Stainless Steel	Yes	-20°C	-	Wall	-

<sup>\*</sup>Ex-proof zone 1 and zone 2.

For more information and technical datasheets please contact your BROEN-LAB sales representative.

-20°C

Outdoor pedestal mounted emergency eye/facewash fountain, trace tape heated and pre-insulated. Delivered with "unbreakable" integral lid, aerated diffusers, flow regulator and self-adhesive sign.

Available in Galvanized or Stainless Steel.

### TECHNICAL SPECIFICATIONS

Min. working pressure: 2 Bar Min. working temperature: -20°C

Water Flow: 20 l/min.

Supply voltage: 230VAC, 50/60Hz



### ORDERING INFORMATION

BROEN-LAB item no.	Material	Ex-proof*	Min. working temperature	Mounting	Heated underground extension pipe
H7 31 001 10 00 321	Galvanized Steel	Yes	-20°C	Floor	Available
H7 31 001 10 00 221	Stainless Steel	Yes	-20°C	Floor	Available

<sup>\*</sup>Ex-proof zone 1 and zone 2.

### **MOBILE & PORTABLE EMERGENCY SHOWERS**

Mobile Emergency Safety Shower is designed as a quick response, first-aid device for use where a plumbed-in unit is not available and must therefore be backed up by a method of further decontamination.

### Mobile Self Contained Emergency Safety Shower



### **Portable Emergency Eye, Face and Body Shower** (Unheated)



No. H7170011000303

### Mobile Self Contained 1200 L Shower and Eyewash incl. Pump

No. H7470031000331

**BROEN-LAB** EMERGENCY SHOWER

mi

### **TECHNICAL SPECIFICATIONS**

Capacity: 114 l

Min. working temperature: -10  $^{\circ}$ C

Water Flow: 40 l/min.

Supply voltage: 230VAC, 50/60Hz

### TECHNICAL SPECIFICATIONS

Capacity: 15 l

Min. working temperature: +1°C

Water Flow: 8 I/min.

Material: Stainless Steel cylinder

### ORDERING INFORMATION

BROEN-LAB item no.	Material	Ex-proof*	Min. working temperature
H7 47 002 10 00 311	Galvanized & SS	No	-10°C
H7 47 002 10 00 321	Galvanized & SS	Yes	-10°C
H7 47 002 10 00 301	Galvanized & SS	No	+1°C (unheated)

<sup>\*</sup>Ex-proof zone 1 and zone 2.



### TECHNICAL SPECIFICATIONS

Capacity: 1200 l

Min. working temperature: -10°C

Water Flow: 75 I/min.

Supply voltage: 230VAC, 50/60Hz

Material:

- Galvanized mild steel trailer

- Red polyethylene tank

A range of heavy duty emergency tank showers, for use when warm water is required, or a constant supply of water cannot be guaranteed.

These units are essential where there is a problem with insufficient water pressure or inadequate flow rate. This improved range incorporates models to suit a wide range of situations.

### Body Shower with Eye Wash Fountain & 350 I tank



### Body Shower with Eye Wash Fountain & 1200 I tank



### **Standard Configuration**

- 350 or 1200 l tank in medium density polythene
- Support legs and frame in galvanized or Stainless Steel
- Rear and side GRP panelling
- Visual level indicator
- Body shower
- · Eye/facewash fountain
- Shower control panic bar
- · Foot panel for body shower

### **TECHNICAL SPECIFICATIONS**

Min. working temperature: -35°C

Water Flow: 75 l/min.

Supply voltage: 230VAC, 50/60Hz

### Tank showers are available in three sizes:

The **350 liter** model supplies a guaranteed flow of warm water for several minutes.

With the mains supply permanently connected, the **750\* liter** model will deliver warm water for 15 minutes, recommended by ANSI Z 358.1 and EN15154-1&2 standards.

The **1200 liter** models will meet the mentioned standards even without an external supply connected.

Available in Galvanized or Stainless Steel

For areas with high outside temperatures models with chillers can be supplied.

Unheated models are also available.

### ORDERING INFORMATION

BROEN-LAB item no.	Material	Ex-proof*	Min. working temperature	Eye wash fountain	Tank volume
H7 61 001 10 01 341	Galvanized Steel	Yes	-35°C	Yes	350 l
H7 61 001 10 01 241	Stainless Steel	Yes	-35°C	Yes	350 l
H7 61 001 10 01 331	Galvanized Steel	No	-35°C	Yes	350 l
H7 61 001 10 01 231	Stainless Steel	No	-35°C	Yes	350 l
H7 61 003 10 01 341	Galvanized Steel	Yes	-35°C	Yes	1200 l
H7 61 003 10 01 241	Stainless Steel	Yes	-35°C	Yes	1200 l
H7 61 003 10 01 331	Galvanized Steel	No	-35°C	Yes	1200 l
H7 61 003 10 012 31	Stainless Steel	No	-35°C	Yes	1200 l

<sup>\*</sup>Ex-proof zone 1 and zone 2.

For more information and technical datasheets please contact your BROEN-LAB sales representative.

<sup>\*750</sup>l tank shower is available on request

Heated underground extension pipe, trace tape heated and pre-insulated.
To be used for connecting to an underground supply.
Available in Galvanized or Stainless Steel.

### TECHNICAL SPECIFICATIONS

Min. working temperature: -35°C Supply voltage: 230VAC, 50/60Hz

### ORDERING INFORMATION

BROEN-LAB item no.	Material	Ex-proof*	Min. working temperature	Mounting
H7 51 001 10 00 321	Galvanized Steel	Yes	-35°C	Underground
H7 51 001 10 00 221	Stainless Steel	Yes	-35°C	Underground

<sup>\*</sup>Ex-proof zone 1 and zone 2.

### **Foot Panel for Operation of Heated Eye Wash**



Weight: Approx. 23 kg

No. 17 803 009 - Stainless Steel version
No. 17 833 009 - Galvanized Steel version

Foot panel or treadle for operation of *heated body shower* is also available - contact your BROEN-LAB sales representative for further information.

### Foot Treadle for Operation of Heated Eye Wash



No. 17 804 009 - Stainless Steel version

For heated eye washes operational down to -20°C

No. 17 805 009 - Stainless Steel version

For heated eye washes operational down to -35°C

For more information and technical datasheets please contact your BROEN-LAB sales representative.

BROEN-LAB EMERGENCY SHOWER SYSTEMS

### **Emergency Showers Are Important**

Thousands of industrial accidents happen every year in spite of precautionary safety measures, accidents that are caused by human exposure to heat (fire, boiling water and oil, steam, etc.), to corrosive agents (fluid acids and caustics) and/or to influence by dust and corrosive aerosols (metal, dirt, fluids).

These accidents cause the injured person acute pain, disablement and may even lead to loss of life.

However, experience shows that the consequences and human suffering following such accidents can be considerably limited if water is applied directly to the affected areas instantly after the exposure.

It is a fact of life that the risk lurks in many places: in the laboratory, in the kitchen, in the warehouse, in the chemical plant, in the workshop or generally speaking in any place where heat, corrosion or dust can cause serious human injury.

For these reasons we strongly recommend the installation of adequate emergency shower systems as a means of limiting these risks.

### **Location and Operation of Emergency Showers**

As a basic rule the shower must be located as close to the workplace as possible so that the general first aid instructions of flushing the affected area immediately after the accident has occurred can be followed.

Please note that if strong acids or caustics are involved immediate action is demanded.

When installing the emergency shower system please also be aware that it must be easy to locate and access at all times for potential users.

The BROEN-LAB showers can be recognized by their red body shower pull rods and red eye shower push handles. In addition the internationally adapted graphic symbols for the clear identification of emergency shower installations are always included in BROEN-LAB solutions.

Further information and guidelines can be found in the EN 15154 part 1 & 2: 2006 - European standard for Emergency safety showers

### Operation

ΑII **BROEN-LAB** emergency designed shower systems are and secure engineered to easy operation. This allows for "one-action" activation of the shower system ensuring immediate and continued flushing without any need for further action.

This easy operation also allows the injured person to use both hands to keep his eyes open for effective flushing, to undress under the shower or just to help exposing the affected area for better and more expedient treatment.

On the eye shower the nozzles are protected by rubber cups to keep the injured person from hurting himself unintentionally. The nozzles are also protected against dust and dirt by the caps which are removed automatically by the initial jet of water.

### **Water Pressure and Temperature**

Emergency showers shall be connected to a cold water system with at least drinking water quality which ensures clean and fresh water. In addition we recommend the water temperature to be kept within the range of 15° to 25° Celsius in order to provide the best immediate effect, reasonable human comfort and at the same time avoid hypothermia.

All eye showers are delivered with FLOWFIX water regulating devices ensuring constant adequate water flow regardless of high water pressure.

The illustrations will inform you of the capacities of the showers.

### Threading

In general the BROEN-LAB Emergency Showers are equipped with threadings as follows:

E.g. G 1/2": In accordance to ISO 228/1.
E.g. NPT 1/2": In accordance to
ANSI-B1.20.1.

### **Maintenance and Information**

When installing emergency shower systems, the relevant safety advisor must provide instructions on operation and maintenance!

As basic information the potential user must be advised of:

- The location of the emergency shower system
- · How to operate it and

once a week.

 How to behave in case of an emergency BROEN-LAB recommends users to check installed emergency shower systems

### IN CASE OF AN EMERGENCY

In the event of accidents caused by human exposure to heat, corrosive agents and/or to influence of dust and corrosive aerosols:

- Start flushing immediately the seconds count.
- Remove contaminated clothing which is not stuck to the skin.
- Continue flushing.
- In case of an eye injury keep the eyelid open when continuing flushing.
- Call for an ambulance. Ask the doctor for advice.
- Continue flushing in the ambulance.
- Keep the flushing in the emergency room and during transport to the special ward.

In case of blistering do inform the doctor and the hospital about the substance in question and about possible countermeasures.

RAPID ACTION CAN BE OF VITAL IMPORTANCE.

### PERSONAL WASH UNIT

The first seconds following an eye injury are often critical to keeping eye injury to a minimum. A personal wash unit may be kept in the immediate vicinity of employees working in a potentially hazardous area.

The main purpose of these units is to supply immediate flushing. With this accomplished, the injured individual should then proceed to a plumbed or selfcontained eyewash and flush the eyes for the required 15-minute period.

REFERENCE: A1 PERSONAL WASH UNIT, ANSI Z358.1-2014


### INTERNATIONAL STANDARDS

### EN-NORM 15154-1&2

(European Committee for Standardization)

### **MAIN REQUIREMENTS**

- The water supplied by the body shower shall be of constant flow rate in accordance with national regulations at a flow pressure to be specified by the manufacturer. The flow pressure shall be specified and measured where the shower is connected to the water system. The body shower shall be capable of delivering this supply for a minimum of 15 minutes. Note: When no national or local regulations apply, a constant flow rate of min. 60 l/min is suitable.
- It shall only be possible to make adjustments with a tool to the direction of spray or the water distribution of a shower head.
   The shower head shall be self-draining between the valve and the outlet. The shower head shall be removable for maintenance but only by use of a tool.
- At a distance 700 mm below the shower head, 50% (+/-10%) of the volume of water delivered shall fall in a circle with a radius of 200mm; the water level in the individual compartments in this circle shall not deviate by more than 30% from the mean value. At this measuring level, the area reached by minimum 95% of the water shall be limited to a circle with a radius of 400 mm. The velocity of the water spray shall be low enough to be non-injurious to the user.
- The shower head shall be designed to be installed so that its lower edge is (2200 +/- 100) mm above the level on which the user stands.
- The valve shall not close automatically once it has been opened. The direction of operating the valve actuator shall be clearly visible and unmistakable.



Please refer to www.BROEN-LAB.com for a 3rd party certificate for EN 15154

### INTERNATIONAL STANDARDS

### ANSI NORM Z 358.1-2014


(American National Standard Institute)

### **MAIN REQUIREMENTS**

- Emergency eyewash and showers must be located in such a way that they can be reached by the user within maximum 10 seconds and available for immediate use (ref.: 4.5.2 + Safety considerations B5).
- Eye/face showers must be able to deliver at least 11.4 l. water per minute continuously for at least 15 minutes (ref.: 6.1.6).
- The handle of pull rod must be located at a maximum height of 173.3 cm above the floor (ref.: 4.2).
- The valve must go from "off" to "on" position within 1 second or less and be able to remain open without further action. (ref.: 4.2).
- The body shower must be capable of delivering flushing fluid at a minimum of 75.7 liters per minute continuously for a minimum of 15 minutes (ref.: 4.1.4).

- The diameter of the water spray pattern must be minimum 50.8 cm measured at 152.4 cm above the floor (ref.: 4.1.4).
- The face of the showerhead shall be located within a range of 208.3 cm to 243.8 cm from the floor (ref.: 4.4.1).
- The eye/face showers flushing fluid flow pattern must be located within a range of 83.8 cm to 134.6 cm above the floor (ref.: 5.4.4).
- The shower must be activated at least once a week (ref.: 4.6.2).
- Tepid water is required in all showers (ref.: 4.5.6 + Safety considerations - B5).

Reference: ANSI Z 358.1-2014



An eye shower is the only allowable exception from the rule that no obstructions are allowed in the "open space".

Please refer to www.BROEN-LAB.com for a 3rd party certificate for ANSI Z358.1-2014

### **BROEN-LAB** EMERGENCY SHOWER SYSTEMS

### SPECIFICATION OF MATERIALS

Shower outlet Brass/PP or PPO Stand pipe

Stainless steel (AISI 304)

or brass

Rubber cups **EPDM Dust caps POM** 

Bowl Stainless steel (AISI 304) **SELF-DRAINING SHOWER HEAD** 

BROEN-LAB has developed a new self-draining emergency shower head in complience with both EN15154-1 and ANSI Z358.1-2014.

BROEN-LAB has 3rd party certificate for EN 15154 and

ANSI Z358.1-2014. The self-draining shower head requires min. 1 bar (14.5 psi) and max. 10 bar (145 psi) dynamic water pressure for

optimal function.



### Available in chemical resistant "BROEN-LAB Polycoat" (polyester powder coating) or "Stanniloy" plating (a highly anti-corrosive treatment based on a tin/nickel alloy).


### **MIN. WORKING PRESSURE** (DYNAMIC/FLOW PRESSURE)

See individual products.

kPa kN/m²		bar	psi
1000	1000	10	145


**MAX. WORKING PRESSURE (STATICAL)** 

### **FLOWFIX**



FLOWFIX is a feature that is built into most of BROEN-LAB's eye baths and hand/ eye showers. Flowfix is a water regulating device, that controls the water flow at a preset rate regardless of water pressure.

### **FEATURES & BENEFITS**




Complies to both EN & ANSI Standards

Chemical resistant PPO plastics

Soft water spread

Available for 1/2", 3/4" and 1" pipes



### CALCULATING THE FLOWRATE AND DYNAMIC PRESSURE DROP

Any calculation can be done using the following formulas:

$$K_v = \ \frac{q_{v \, [m^3/h]}}{\sqrt{\Delta P_{[bar]}}} \qquad \text{or} \quad q_v = K_v \times \sqrt{\Delta P} \qquad \text{or} \quad \Delta P = \left(\frac{q_v}{K_v}\right)^2$$

### Where

 $K_v = \text{flowrate } [m^3/h] \text{ or } [l/\text{sec}] \text{ at 1 bar pressure drop}$ 

q<sub>v</sub> = flowrate [m<sup>3</sup>/h] or [l/sec]

 $\Delta P$ = pressure drop [bar]

1) Please note that eye showers with built-in FLOWFIX have a fixed flowrate

### **Examples:**

What pressure drop is needed in order to achieve 75 l/min = 1.25 l/sec through a body shower with 1/2" inlet?

$$\Delta P = \left(\frac{q_v}{K_v}\right)^2 = \left(\frac{1.25 \text{ [l/sec]}}{1.08 \text{ [l/sec]}}\right)^2 = 1.34 \text{ [bar]}$$

What flowrate can be expected through a free standing body shower with top inlet at a pressure drop of 1.5 bar?

$$q_v = K_v \times \sqrt{\Delta P} = 1.10 \text{ [l/sec] } x\sqrt{1.5} \text{ [bar]} = 1.35 \text{ l/sec} = 80.8 \text{ l/min.}$$

### **SYMBOLS:**



Eye shower



**Body shower** 

### **K<sub>v</sub> VALUES**

Type o		Eye shower without FLOWFIX <sup>1)</sup>	Body shower 1/2" pipe or inlet	Body shower 3/4" pipe or inlet	Body shower 1" pipe or inlet	Free standing- body shower, top inlet	Free standing- body shower, bottom inlet	Thermostatic mixing valve 17 140.009	Thermostatic mixing valve 17 459.009
K <sub>v</sub> -	m³/h	0.75	3.89	4.3	5.31	3.97	2.99	1.25	5.20
value	l/sec	0.21	1.08	1.19	1.48	1.10	0.83	0.35	1.44



For disabled people as well

NOTES			

NOTES	
	-
	<u></u>
	7

BROEN-LAB develops, manufactures and sells Laboratory Fittings, Emergency Shower Systems and Eye Wash Systems into a broad spectrum of laboratories and industrial locations; sectors include Pharmaceutical, Food&Beverage and Academia. Our expertise and product quality ensure optimal solutions compliant to all relevant international norms setting new standards in risk mitigation in modern research and development facilities.

BROEN-LAB offer solutions that ensure the functionality, compliance, hygiene, durability and safety are of the highest priority. For over 50 years our products have been integrated into a wide variety of workplaces, laboratories, hospitals and industrial locations, with features that are assessed and designed to be flexible, durable and compliant offering a broad range of options tailored to each location. This inherent design integrity provides confidence in the solution adopted wherever in the world these are applied.

BROEN-LAB is a collaboration partner in the design and layout of your next laboratory and it's integrated safety solutions.

### **HEAD OFFICE**

### **BROEN-LAB A/S**

Drejervænget 2 • DK-5610 Assens Tel. +45 6376 6376 • Fax +45 6471 2476 E-mail: lab@broen-lab.com



### **BROEN-LAB A/S ISO 9001 certification**

In September 1991 BROEN-LAB A/S was certified according to ISO 9001 as one of the first Danish companies. The certification was carried out by Bureau Veritas, Denmark, for our Danish site in Assens. The quality management system of BROEN-LAB A/S now complies with detailed specifications laid down by the internationally acknowledged EN ISO 9001:2008. This certification will further contribute to reputation for quality and reliability of BROEN-LAB A/S.

### **INTERNATIONAL OFFICES**

### **BROEN-LAB Ltd.**

Unit 9 Navigation Point, Golds Hill Way Tipton • West Midlands DY4 OPY • U.K. Tel. +44 121 522 4515 • Fax +44 121 522 4535 E-mail: lab@broen-lab.com

### **BROEN-LAB Singapore Pte. Ltd.**

10 Bukit Batok Crescent #10-02 • The Spire • Singapore 658079 Tel: +65 6298 0662 • Fax: +65 6298 0468 E-mail: lab@broen-lab.com

### **BROEN-LAB GmbH**

Tel. +45 6376 6376 • Fax +45 6471 2476 E-mail: lab@broen-lab.com

### **BROEN-LAB Sverige AB**

Tel. +45 6376 6376 • Fax +45 6471 2476 E-mail: lab@broen-lab.com

www.BROEN-LAB.com